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## ENGLISH

# IMPORTANT SAFETY NOTES

READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO ALLOW FOR CORRECT USE OF THIS PRODUCT.

PLEASE FIND THE COMPLETE MANUAL UNDER THE FOLLOWING LINK:  
[www.runtastic.com/products](http://www.runtastic.com/products)

### WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

### CAUTION

CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

### NOTICE


NOTICE addresses practices not related to personal injury, such as product and/or property damage.

Observe the instructions for use.

Type B applied part.

	<p>Disposal in accordance with EC Directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment).</p>
	<p>Means that batteries should not be disposed of with household waste!</p>
	<p>Operating temperature</p>
	<p>Operating humidity</p>
	<p>Protect from moisture</p>
	<p>The CE labeling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products. Distributor (Importer): Beurer GmbH, Söflinger Str. 218, 89077 Ulm, Germany.</p>

## CLASSIFICATION

- a. Internally powered equipment;
- b. Type B applied part;
- c. IPX0;
- d. Not Category AP/APG equipment;
- e. Mode of operation: continuous operation.

## **FCC Compliance Statement (For US)**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: **(1)** this device may not cause harmful interference, and **(2)** this device must accept any interference received, including interference that may cause undesired operation.

### **NOTE:**

This equipment has been tested and found to comply with the limits of a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void user's authority to operate the equipment.

## IC Statement (For Canada):

**EN:** This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: **(1)** this device may not cause interference, and **(2)** this device must accept any interference, including interference that may cause undesired operation of the device.

This Class [B] digital apparatus complies with Canadian ICES-003.

Operating temperature: +41°F to +104°F (+5°C to +40°C)

RH ≤ 95%

Operating humidity: 20 to 95% RH

## **BATTERY**

### **WARNING**

- Always keep batteries out of the reach of children. Batteries can be life threatening if swallowed. Therefore, keep the batteries and scale out of the reach of small children. If a battery is swallowed, medical assistance should be obtained immediately.
- Remove a drained battery from the scale promptly, otherwise it can leak. In this case avoid contact with skin, eyes and mucous membranes. Rinse any areas affected by battery acid immediately with plenty of clear water, and seek medical assistance immediately.

### **CAUTION**

- Batteries may not be charged or reactivated by other means, taken apart, thrown into fire or short circuited.
- If necessary, clean the battery contacts and also the counter contacts in the device.
- Danger of explosion if the battery is replaced incorrectly.
- Batteries are to be inserted with correct polarity.
- Never mix old and new batteries (replace all batteries at the same time).
- Never mix Alkaline, standard or re-chargeable batteries.

- Use only batteries of the same or equivalent types as recommended.
- Do not recharge non-rechargeable batteries.
- Never expose the batteries to excessive heat such as sunshine, fire or similar! Non-adherence to these instructions can lead to damage and sometimes even explosion of the batteries.

## NOTICE

- Clean the batteries' contacts and also those of the product prior to battery installation.
- The supply terminals are not to be short-circuited.
- Always remove exhausted or dead batteries.
- Remove the battery from the scale when not being used for long periods of time.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable local regulation. Do not dispose of batteries with normal household waste.

If you wish to dispose of your scale you must comply with the applicable regulations. Products that are designated with the symbol below must not be disposed of as normal household waste. Details are provided by the community collection point.

## **WARNING**

- The scale should not be used by individuals with medical implants (i.e. pacemakers).  
The function of the implanted device may be impaired.
- The scale is not to be used during pregnancy.
- Not intended for use by children under 10.
- Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- You can find further important safety information in the complete manual which can be found under the following link:  
[www.runtastic.com/products](http://www.runtastic.com/products)

# DECLARATION OF CONFORMITY

EMC Directive 2004/108/EC, R&TTE Directive 1999/5/EC, RoHS Directive 2011/65/EU:

Short text of the conformity declaration: Krippel-Watches hereby declares that the Bluetooth scale (RUNSCA1) complies with the basic requirements and the other relevant provisions of the EMC Directive 2004/108/EC, RoHS Directive 2011/65/EU and the remaining relevant conditions of the R&TTE Directive 1999/5/EC.

You can access the full text of the conformity declaration on our homepage:  
[www.runtastic.com/products](http://www.runtastic.com/products)

# EMC COMPLIANCE DECLARATION

## ELECTROMAGNETIC COMPATIBILITY

### INFORMATION

Table 1

**For all ME EQUIPMENT and  
ME SYSTEMS**

Guidance and manufacture’s declaration – electromagnetic emissions		
The LIBRA is intended for use in the electromagnetic environment specified below. The customer or the user of the LIBRA should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The LIBRA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	

Table 2

**For all ME EQUIPMENT and  
ME SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic immunity			
The LIBRA is intended for use in the electromagnetic environment specified below. The customer or the user of the LIBRA should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3

**For ME EQUIPMENT and ME SYSTEMS  
that are not LIFE-SUPPORTING**

Guidance and manufacturer’s declaration – electromagnetic immunity			
The LIBRA is intended for use in the electromagnetic environment specified below. The customer or the user of the LIBRA should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	[E <sub>1</sub> ] V/m 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the LIBRA, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W)</p>

			according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
--	--	--	--

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LIBRA is used exceeds the applicable RF compliance level above, the LIBRA should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the LIBRA.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Table 4

**For ME EQUIPMENT and ME SYSTEMS  
that are not LIFE-SUPPORTING**

**Recommended separation distances between portable and mobile RF communications equipment and the LIBRA**

The LIBRA is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the LIBRA can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LIBRA as recommended below, according to the maximum output power of the communications equipment.

Rated maximum outputpower of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.117	0.233
0.1	0.368	0.737
1	1.167	2.333
10	3.689	7.37
100	11.667	23.333

---

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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